



463PR

G-Biosciences ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ [technical@GBiosciences.com](mailto:technical@GBiosciences.com)

A Geno Technology, Inc. (USA) brand name

# Yeast PE LB™

## Yeast Protein Extraction Lysis Buffer

(Cat. # 786-178, 786-179)



think proteins! think G-Biosciences [www.GBiosciences.com](http://www.GBiosciences.com)

INTRODUCTION ..... 3

APPLICATIONS ..... 3

COMPATIBILITY ..... 3

ITEM(S) SUPPLIED ..... 3

STORAGE CONDITIONS ..... 3

ADDITIONAL ITEMS NEEDED ..... 4

PREPARATION BEFORE USE ..... 4

PROTOCOL ..... 4

    SPHEROPLAST LYSIS ..... 5

RELATED PRODUCTS ..... 5

## INTRODUCTION

Yeast PE LB™ is useful for extraction of soluble proteins from yeast cells. Yeast PE LB™ is a proprietary improvement on the Zymolyase® based spheroplast preparation and extraction of soluble proteins from yeast cells. This kit is provided with an optional protocol to make spheroplast and remove lytic enzyme Zymolyase®, prior to lysis and extraction of yeast proteins. Yeast PE LB™ is based on organic buffering agents that utilize a mild non-ionic detergent and a proprietary combination of various salts and agents to enhance extraction and stability of proteins. A ready-to-use Zymolyase® preparation is also provided. Depending on application, additional agents such as reducing agents, chelating agents, and protease inhibitors may be added into Yeast- PE LB™ (see Related Products for protease inhibitor Protease Arrest™). The proprietary combination of this reagent provides a simple and versatile method of yeast protein extraction. Yeast PE LB™ eliminates the need for laborious glass bead lysis of yeast cells. The kit is good for 100 preps of 50µl yeast cell pellet.

## APPLICATIONS

Preparation of yeast spheroplast and extraction of yeast proteins. This kit is suitable for processing approximately 10ml yeast cell pellet suspension, either single or multiple smaller preps.

## COMPATIBILITY

Yeast PE LB™ is compatible with any downstream application including running various chromatography procedures and gel electrophoresis applications. Yeast PE LB™ is also compatible for protein estimation with NI™ protein assay.

## ITEM(S) SUPPLIED

Description	Cat. # 786-178	Cat. # 786-179
Yeast PE LB™ Buffer	100ml	500ml
Yeast Suspension Buffer	15ml	-
Longlife™ Zymolyase® (1500U/ml)	2 x 0.5 ml	-

## STORAGE CONDITIONS

The kit is shipped at ambient temperature. Upon arrival store the kit components at 4°C except Longlife™ Zymolyase® at -20°C. Stable for 1 year when stored and used as recommended.

## ADDITIONAL ITEMS NEEDED

- Centrifuge
- Test tubes
- Incubator,
- DTT, EDTA and  $\beta$ -mercaptoethanol

Additional volume of the Yeast PE LB™ Buffer may be purchased separately for downstream applications such as chromatography, dialysis, etc.

## PREPARATION BEFORE USE

*Depending on applications, DTT and EDTA may be added. Prepare an appropriate volume of the Yeast PE LB™ for use by adding DTT and EDTA both to a final concentration of 5mM. If the presence of a divalent metal ion is necessary for any application, do not add EDTA; instead add an appropriate divalent salt to a final concentration of 5mM.*

*Protease Inhibition: If the inhibition of protease activity is required, add a cocktail of protease inhibitors to prevent protease activities during extraction procedure We recommend our ProteaseARREST™ Protease Inhibitor Cocktail (Cat. # 786-108).*

## PROTOCOL

1. Pellet Yeast cells (culture OD<sub>600</sub> 1.5-2.0) by centrifugation at 5-10,000x g for 10 minutes. Suspend the cell pellet in an equal volume of the Yeast Suspension Buffer. Add 1 $\mu$ l of  $\beta$ -mercaptoethanol per 100 $\mu$ l Yeast suspension.
2. Vortex for 1 minute or until the cell suspension is homogeneous. Incubate the cell suspension for 5 minutes at 4°C. Vortex it again to suspend the cells.
3. Flick the vial containing Longlife™ Zymolyase® to mix the solution. Add 10 $\mu$ l Longlife™ Zymolyase® for each 100 $\mu$ l cell suspension. Gently mix the content.
4. Incubate the suspension at 37°C for 30-60 minutes. Lysis can be monitored by taking 25 $\mu$ l suspension, mixing with 1ml Yeast PE LB™ Buffer and reading optical density at 800nm.
5. At the end of incubation, centrifuge the suspension at 10,000x g for 5 minutes. Remove and discard the supernatant carefully, leaving the spheroplast pellet in the tube.

***Optional Washing:*** Add 10 volumes (~600 $\mu$ l) of the Yeast Suspension Buffer to the spheroplast pellet. Re-suspend the spheroplast by gently tapping the tube. Centrifuge again as above and discard the supernatant.

## **Spheroplast Lysis**

1. Suspend the yeast spheroplast pellet in 3 volumes of the Yeast PE LB™ Buffer. Pipette the suspension up and down a few times. Vortex periodically and incubate on ice for 30 minutes. Incubating the cells for 3 minutes at 37°C or a brief sonication step may further facilitate the lysis. Sonication is necessary for shearing genomic DNA.

**NOTE:** For better lysis, increase the Yeast PE LB™ to spheroplast ratio, because higher the ration better would be the lysis.



2. Centrifuge at 20,000x g for 30 minutes at 4°C. Collect clear lysate. The lysate is now ready for purification of protein, other applications, or further analysis.

**NOTE-** Additional volume of Yeast PE LB™ can be purchased separately for downstream applications e.g. chromatography and dialysis, etc.

Zymolyase® is a registered trademark of Kirin Brewery Co. Ltd.

## **RELATED PRODUCTS**

Download our Sample Preparation and Protease & Phosphatase Inhibitors, Enzymes & Assays Handbooks



<http://info.gbiosciences.com/complete-protein-sample-preparation-handbook/>  
<http://info.gbiosciences.com/protease-phosphatase-inhibitors-enzymes-assay-handbook/>

For other related products, visit our website at [www.GBiosciences.com](http://www.GBiosciences.com) or contact us.

Last saved: 9/11/2012 CMH

*This page is intentionally left blank*

*This page is intentionally left blank*



[www.GBiosciences.com](http://www.GBiosciences.com)