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文章來源：The benefits of using capacitance as a direct measurement of cell density and health in bioprocessing (ABER)

## References

- [1.Kell D B. DETERMINATION OF BIOMASS World Patent WO1988002114 \(1988\)](#)
- [2.Butler M, et al. Cell Viability in Bioprocesses Making a Case for Reevaluation. BioProcess International. 2018: 17 \(11–12\): 42-49](#)
- [3.Carvell J P, Dowd J E. On-Line Measurement and Control of Viable Cell Density in Cell Culture Manufacturing Processes Using Radio-Frequency Impedance. Cytotechnol. 2006: 50 \(1–3\): 35–48.](#)
- [4.Metze, S et al. Monitoring online biomass with a capacitance sensor during scale-up of industrially relevant CHO cell culture fed-batch processes in single-use bioreactors. Bioprocess and biosystems engineering. 2020: 43 \(2\): 193-205.](#)
- [5.Carvell J, et al. Monitoring Live Biomass in Disposable Bioreactors. BioProcess International 2016: 14\(3\)s: 40-48.](#)
- [6.Madsen B, et al. Measuring Cell Density in HyPerforma S.U.B.s with ABER Futura neotf. Bioprocess International 2021 special report:](#)
- [7.Fernandes J, et al. Development of capacitance tools: at-line method for assessing biomass of mammalian cell culture and fixed cell calibration standard. Biotechnol J. 2018: 14 \(4\): e1800283.](#)
- [8.Braasch K, et al. The changing dielectric properties of CHO cells can be used to determine early apoptotic events in a bioprocess. Biotechnol Bioeng. 2013:110 \(11\):2902-2914.](#)
- [9.Ma F et al. Real-time monitoring and control of CHO cell apoptosis by in situ multifrequency scanning dielectric spectroscopy. Process Biochemistry, 2019: 80: 138-145.](#)
- [10.Moore, B. et al. Case study: The characterization and implementation of dielectric spectroscopy \(biocapitance\) for process control in a commercial GMP CHO manufacturing process. Biotechnol. Prog. 2019:35\(3\): e2782.](#)
- [11.Bro, C. Kwiatkowski, C and Tolstrup, A. Use of biocapitance probes for optimized process control at large-scale manufacturing in Cell Culture Engineering XVI, A. Robinson, PhD, Tulane University R. Venkat, PhD, MedImmune E. Schaefer, ScD, J&J Janssen Eds, ECI SymposiumSeries, \(2018\).](#)
- [12.Zhang, A., et al. Advanced process monitoring and feedback control to enhance cell culture process production and robustness. Biotechnol. Bioeng. 2015: 112: 2495-2504.](#)
- [13.Morris, C et al. Single in-line biomass probe detects CHO cell growth by capacitance and bacterial contamination by conductivity in bioreactor. Biotechnol J. 2021 S30: e2100126.](#)